

## EEBCSTAT SUPPORT PANEL WITH INSTRUCTIONS



SKU: 149633F



**Wall mounting support panel with English instructions** for EEBCSTAT test stations.

Includes 2 cables for the connection to tester and footwear plate.

**Text on the panel:** "Test Station. Wrist-Strap Test (Wearing Strap). Connect Wrist-Strap Lead To Reference Point. Press Test Button. Footwear Test (Wearing Heel Strap Or Toe Strap) Connect Foot Plate To Reference Point. Stand On Footplate. Press Test Button. Reference Point".

**Standard supply:** EEBCSTAT support panel, Connection cable approx. 15 cm (0 mega-Ohm), Connection cable approx. 182 cm (0 mega-Ohm).